

Previously known as: MICROGARD® 2500 Overboots SOCO - Model 409

Durable and breathable overboots, designed with microporous laminate technology.

- **Certified performance:** Developed specifically for police forensic scene of crime officers, AlphaTec™ 2500 STANDARD Overboots SOCO Bound - Model 409 complies with EN 14126 standards and affords superb protection from biological agents, including blood, bodily fluids and blood-borne pathogens
- **Advanced defenses:** Owing to its unique microporous polypropylene laminate fabric, this durable yet comfortable personal protection footwear guards against harmful chemical sprays and infective agents
- **Reduced contamination risks:** In addition, these chemical-resistant overboots are particularly suited for critical areas, as they are low linting and do not pose any contamination risks



Industries

- Defense


Applications

- Transferring liquids and solids
- Spill or leakage cleanup
- Unexpected leakages spills or other releases
- Plant disinfection and sanitization
- Maintenance
- Emergency Response
- Contaminated Environments


KEY FEATURES & BENEFITS

- ◊ **EN-certified protection:** Suitable for even the harshest environments
- ◊ **Unique material:** Biological and chemical resistance assured
- ◊ **Low linting:** Lower contamination risks, increased product protection

Material Diagram



Seam Diagram



Performance Standards



PARTIAL BODY
PROTECTION ONLY



TYPE PB [6]



0598

EN ISO 13688:2013



EN 14126



Product Information

Product Information	
Available Sizes	One Size
Color	White
Packaging Overview	20 pairs every compressed bag
Product Material	Microporous Polypropylene Laminate non-woven
Product Reference	2500 Overboots SOCO - Model 409
Standards Overview	CE CATEGORY III (Type PB (6)-B, EN 14126)
Seam Type	Bound
Shelf Life	5 years
Style Number	W25B409

Suit components & accessories

Suit Components and Accessories	
Design Features	SOCO Overboots, tie fastening, blue binding to seams, reinforced Surestep non-slip soles, adjustable shoe tie.

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region
Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00

North America Region
Ansell Healthcare Products LLC
T: +1 800 800 0444

Australia
Ansell Limited
T: +61 1800 337 041

Asia Pacific Region
Ansell Global Trading Center
T: +603 8310 6688

Latin America & Caribbean Region
Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

