



## Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

# Sterile polychloroprene cleanroom gloves, offering users chemical resistance and tactility

- **Specialized protection:** BioClean™ Ultimate BUPS cleanroom gloves are tested against ASTM D6978 and EN standards, guarding against chemotherapy drugs and various chemicals
- **Assured component integrity:** These latex-free gloves' electrostatic-dissipative (ESD) properties ensure greater component protection against electrostatic discharge
- **Reduced contamination risks:** Their non-particulating EasyTear™ packaging is designed to minimize contamination risks
- **Enhanced features:** These powder-free cleanroom gloves are manufactured with a beaded cuff, enhancing their strength and facilitating double-donning
- **Heightened tactility and comfort:** BioClean™ Ultimate BUPS polychloroprene gloves are also comfortable, with good tactility levels for engaging in precision work



## KEY FEATURES & BENEFITS

- ◊ **ASTM and EN standards:** Chemical and chemotherapy drug resistance
- ◊ **ESD properties:** Component integrity protection
- ◊ **Non-particulating EasyTear packaging:** Fewer contamination risks

## Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Sterile Processing Department
- Lab and R&D



[WWW.ANSELL.COM](http://WWW.ANSELL.COM)

**Ansell**



# Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

## TECHNICAL DATA SHEET

|                                        | Product Information                                                                                                                                 |
|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                               | Neoprene (Polychloroprene)                                                                                                                          |
| Color                                  | Natural                                                                                                                                             |
| Shape                                  | Hand Specific                                                                                                                                       |
| Cuff                                   | Beaded                                                                                                                                              |
| Manufacturing/QMS Audit Standards      | ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms                                  |
| Regulatory/Standard Compliance         | ASTM D6978, CE 0598, EN ISO 21420:2020, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Category III                              |
| Packaging                              | 1 pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs) |
| Storage                                | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.                                                                 |
| Available sizes                        | 6, 6.5, 7, 7.5, 8, 8.5, 9                                                                                                                           |
| External Glove Surface                 | Textured Fingers and Palm                                                                                                                           |
| Internal Glove Surface                 | Polymer Coated                                                                                                                                      |
| Sterilization Method                   | GAMMA irradiation (25 kGy)                                                                                                                          |
| Sterilization Minimum Dose             | 25kGy                                                                                                                                               |
| Sterility Assurance Level              | 10 <sup>-6</sup>                                                                                                                                    |
| Cleanroom Class                        | Class 10/ISO Class 4, EU GMP Grade A/B and other sterile cleanrooms                                                                                 |
| Shelf Life                             | 5 Years                                                                                                                                             |
| Tested for use with Chemotherapy Drugs | Yes, in accordance with ASTM D6978                                                                                                                  |
| Protein Level                          | Not Applicable: contains no natural rubber latex                                                                                                    |
| Anti-static                            | Yes                                                                                                                                                 |



# Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

| Physical Properties                             |                              |        |        |        |       |         |         | Testing Method  |
|-------------------------------------------------|------------------------------|--------|--------|--------|-------|---------|---------|-----------------|
| Sizes                                           | 6                            | 6.5    | 7      | 7.5    | 8     | 8.5     | 9       | EN 420          |
| Typical Length (mm/in)                          | 300 / 12                     |        |        |        |       |         |         |                 |
| Palm Width (mm/in)                              | 77/3                         | 83/3.2 | 89/3.5 | 95/3.7 | 102/4 | 108/4.2 | 115/4.5 |                 |
| Freedom from Holes                              | 0.65 AQL Performance Level 3 |        |        |        |       |         |         | EN 374-2        |
| Typical Particle Count<br>≥0.5µm (counts / cm²) | <1300                        |        |        |        |       |         |         | IEST-RP-CC005.4 |
| Force at Break (N)<br>Before Aging              | ≥ 9 N                        |        |        |        |       |         |         | EN 455-2        |
| Force at Break (N)<br>During Aging              | ≥ 9 N                        |        |        |        |       |         |         | EN 455-2        |

## IONIC CONTENT

| Concentration in µg/cm <sup>2</sup> | Typical      | Concentration in µg/cm <sup>2</sup> | Typical      |
|-------------------------------------|--------------|-------------------------------------|--------------|
| Ammonium                            | 0.005        | Nitrate                             | 0.795        |
| Bromide                             | Not Detected | Nitrite                             | Not Detected |
| Calcium                             | 0.582        | Phosphate                           | Not Detected |
| Chloride                            | 0.494        | Potassium                           | 0.378        |
| Fluoride                            | Not Detected | Sodium                              | 0.123        |
| Lithium                             | Not Detected | Sulphate                            | 0.041        |
| Magnesium                           | 0.025        | Zinc                                | Not Tested   |

## ORDERING INFORMATION

| SIZE        | 6       | 6.5     | 7       | 7.5     | 8       | 8.5     | 9       |
|-------------|---------|---------|---------|---------|---------|---------|---------|
| REORDER NO. | BUPS-60 | BUPS-65 | BUPS-70 | BUPS-75 | BUPS-80 | BUPS-85 | BUPS-90 |

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

