

Long cuff cleanroom gloves, providing extra protection coverage against chemicals and chemotherapy drugs

- **Chemical protection:** BioClean™ N-Plus BNPS are resistant to a range of chemicals and chemotherapy drugs (PPE Category III, EN ISO Type B, ASTM D6978)
- **Clean processed:** The gloves undergo deionized washing to further reduce surface particulates, ionic extractables and non-volatile residues
- **Reduce Contamination Risk:** The gloves are packed in non-particulating packaging reducing potential contaminant hazards in the cleanroom environment



KEY FEATURES & BENEFITS

- ◊ **Extended Long Cuff:** Provide increased coverage for excellent protection of the forearms
- ◊ **Chemical protection:** Resistant to a range of chemicals and chemotherapy drugs
- ◊ **Sterility Assurance:** Gamma sterilized with sterility assurance level (SAL) of 10^{-6}

Industries

- Biotechnology Manufacturing
- Controlled and Critical Environments
- Laboratory Research and Development
- Medical Device Manufacturing
- Pharmaceutical Manufacturing
- Production and Manufacturing
- Sterile Processing Department

Recommended For

- Transferring liquids and solids
- Chemical Splash
- Sample taking and processing
- Biotechnologies Blending operations
- Aseptic filling and labelling



N-Plus BNPS

Sterile Disposable Nitrile Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Nitrile
Color	White
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	ASTM D6978, CE 0493, EN 421:2010, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Contact, Category III, UKCA
Packaging	1 pair per inner wrap; 1 inner wrap per pouch; 10 pouches per outer bag; 20 outer bags per lined carton (200 pairs) More sustainable packaging: These sterile nitrile gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status of PE as this plastic material may not be considered suitable for recycling in your country.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9, 10
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Chlorinated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO Class 4, EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	5 Years
Tested for use with Chemotherapy Drugs	Yes
Protein Level	Not Applicable: contains no natural rubber latex
Anti-static	Yes



N-Plus BNPS

Sterile Disposable Nitrile Cleanroom Glove

Physical Properties									Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	10	EN ISO 21420
Typical Length (mm/in)	400 / 15.75								
Palm Width (mm/in)	80/3.1	85/3.3	93/3.6	98/3.8	106/4.1	110/4.3	115/4.5	120/4.7	
Freedom from Holes	1.0 AQL								EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<1800								IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.17/6.69								EN ISO 21420
Target Single Wall Finger Thickness (mm/mil)	0.21/8.26								EN ISO 21420
Target Single Wall Cuff Thickness (mm/mil)	0.10/3.93								EN ISO 21420
Force at Break (N) During Aging	≥ 9 N								EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	Available On File	Nitrate	<0.30
Bromide	<0.01	Nitrite	Available On File
Calcium	<0.50	Phosphate	<0.01
Chloride	<0.45	Potassium	<0.20
Fluoride	<0.01	Sodium	<0.10
Lithium	<0.01	Sulphate	<0.25
Magnesium	<0.01	Zinc	< 0.18

ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9	10
REORDER NO.	BNPS60	BNPS65	BNPS70	BNPS75	BNPS80	BNPS85	BNPS90	BNPS100

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

