

# ENCORE™ Non-Latex PI Hybrid

S U R G I C A L   G L O V E S

## A hybrid surgical glove ideally balancing polyisoprene (PI) comfort and durability of neoprene

- **Optimal comfort and durability:** These surgical gloves' feature Ansell's patented hybrid blend technology enabling wearers to enjoy the best of two synthetic glove materials: the latex-like comfort of PI and the superior durability of neoprene
- **Heightened allergy prevention:** Being a synthetic glove, ENCORE® Non-Latex PI Hybrid eliminates the risk of a latex allergy while a range of allergens\* known to cause allergic contact dermatitis are eliminated
- **Enhanced safety:** When used as an outer glove for double gloving, this semi-transparent glove allows for fast and easy breach detection
- **Strengthened grip:** These hybrid surgical gloves feature a smooth surface with a micro-textured finish for secure instrument handling
- **Chemo-tested:** These sterile gloves comply with ASTM D6978 standards, making them suitable for use with various chemotherapy drugs\*\*

\* Free from natural rubber latex proteins and do not contain Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and cetylpyridinium chloride (CPC)

\*\* Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations


### Recommended For

- Easy identification of a glove breach
- Protection from Type I latex allergy in HCW's or patients
- Use as an outer glove when double gloving

### FEATURED TECHNOLOGIES

 **A.R.T.**  
Pinhole Barrier Technology

 **DERMASHIELD**  
Donning Aid Technology



## ENCORE® The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

# ENCORE™ Non-Latex PI Hybrid

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                                                                   |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                                | PI (Polyisoprene), Neoprene (Polychloroprene)                                                                                                                     |
| Powder Content                          | Powder-Free                                                                                                                                                       |
| Color                                   | Natural                                                                                                                                                           |
| Thickness                               | Standard                                                                                                                                                          |
| Cuff Style                              | Beaded with SUREFIT™ Technology                                                                                                                                   |
| External Glove Surface                  | Smooth with Micro-Textured Finish                                                                                                                                 |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD™ Technology                                                                                                                       |
| Allergy Prevention                      | Latex (Type I)                                                                                                                                                    |
| Grip Level                              | Moderate                                                                                                                                                          |
| Double Gloving Recommendation           | Outer glove                                                                                                                                                       |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                                                |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                                                                          |
| Sterile                                 | Yes                                                                                                                                                               |
| Product Standards                       | ASTM D6978, ASTM F739, CE 0493 (PPE) + CE 2797 (MD) , EN 16523-1, EN 374-2 and -4, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                                                                                                          |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                                                              |
| Shelf Life                              | 3 years                                                                                                                                                           |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                          |
| Vulcanization Chemical Accelerators     | Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC), Sodium di-n-butylthiocarbamate (SDBC), Xanthogen Polysulphides                            |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 340105055 | 340105060 | 340105065 | 340105070 | 340105075 | 340105080 | 340105085 | 340105090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Ansell Commercial Mexico**  
 Tel: +52 442 296 2050  
 Email: latam.medical@ansell.com  
**Ansell Brasil Ltda**  
 +55 (11) 3356-3100  
 Email: latam.medical@ansell.com  
**Ansell Colombia SAS**  
 Tel: +57 1 288 3247  
 Email: latam.medical@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com