

# GAMMEX™ Non-Latex

## SURGICAL GLOVES

### Neoprene, chemical accelerator-free surgical gloves offer advanced allergy protection

- **Advanced allergy protection:** GAMMEX® Non-Latex neoprene surgical gloves are synthetic and eliminate risks of Type I latex allergy. These sterile gloves are made 100% free of chemical accelerators and minimize any risk of sensitization to allergens leading to allergic contact dermatitis
- **Proven durability:** Strong neoprene glove film deliver the added hand protection needed for the most rigorous procedures
- **Confident grip:** These sterile gloves feature an external micro-textured finish for firmer, more confident instrument handling
- **Effortless donning:** Ansell's donning aid technology includes an internal polymer coating that facilitate easy dry or damp donning
- **Chemo-tested:** In addition, these long-cuff gloves are tested against ASTM D6978 and recommended for the safe handling of chemotherapy drugs\*



\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.

### Recommended For

- Damp donning
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Tested for use with chemotherapy drugs

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX™ Non-Latex

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                                                                                            |
|-----------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                                | Neoprene (Polychloroprene)                                                                                                                                                                 |
| Powder Content                          | Powder-Free                                                                                                                                                                                |
| Color                                   | Green                                                                                                                                                                                      |
| Thickness                               | Standard                                                                                                                                                                                   |
| Cuff Style                              | Straight with SUREFIT™ Technology                                                                                                                                                          |
| External Glove Surface                  | Smooth with Micro-Textured Finish                                                                                                                                                          |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD™ Technology                                                                                                                                                |
| Allergy Prevention                      | Latex (Type I), Chemical (Type IV)                                                                                                                                                         |
| Grip Level                              | Moderate                                                                                                                                                                                   |
| Double Gloving Recommendation           | Outer glove or underglove                                                                                                                                                                  |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                                                                         |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                                                                                                   |
| Sterile                                 | Yes                                                                                                                                                                                        |
| Product Standards                       | AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP                                                                                                                       |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                                                                                       |
| Shelf Life                              | 3 years                                                                                                                                                                                    |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                                                   |
| Vulcanization Chemical Accelerators     | None                                                                                                                                                                                       |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 340006055 | 340006060 | 340006065 | 340006070 | 340006075 | 340006080 | 340006085 | 340006090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Ansell Commercial Mexico**  
 Tel: +52 442 296 2050  
 Email: latam.medical@ansell.com  
**Ansell Brasil Ltda**  
 +55 (11) 3356-3100  
 Email: latam.medical@ansell.com  
**Ansell Colombia SAS**  
 Tel: +57 1 288 3247  
 Email: latam.medical@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com