

# GAMMEX™ PI Radiation Attenuation

SPECIALTY GLOVES

## Polyisoprene sterile gloves, protecting hands from risks of scattered radiation exposure

- ◊ **Enhanced scattered radiation protection:** Suitable for use with C-Arms, these lead-free gloves reduce skin-dose exposure to scattered radiation by 63.0%, compared to just 1%\*\* with standard surgical gloves
- ◊ **Assured dexterity and tactile sensitivity:** Ansell's proprietary, soft, malleable polyisoprene formulation ensures the dexterity, tactile sensitivity and comfort needed for image-guided interventional procedures (fluoroscopy)
- ◊ **Reduced allergy risks:** Being a non-latex surgical glove, GAMMEX® PI Radiation Attenuation eliminates the risk of a latex allergy
- ◊ **Strengthened grip:** These black specialty gloves also feature a textured finish for good grip control




\*\*Wagner LK, Mulhern OR. Radiation-attenuating surgical gloves: effects of scatter and secondary electron production. Radiology.1996;200(1):45-48

CAUTION: These radiation attenuation gloves are designed to protect hands from scattered radiation exposure originating from the X-ray beam during fluoroscopic procedures. They are not intended for use in direct or primary X-ray beams.

### Recommended For

- Protecting hands from image guided interventional procedures (fluoroscopy)
- Protection from Type I latex allergy in HCW's or patients



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



# GAMMEX™ PI Radiation Attenuation

SPECIALTY GLOVES

## PRODUCT INFORMATION

Material	Polyisoprene with Tungsten
Powder Content	Powder-Free
Color	Black
Thickness	Standard
Cuff Style	Straight
External Glove Surface	Textured Fingertips
Internal Glove Surface	Chlorinated
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM F2547, EU: CE Marked to Medical Device Regulation (EU) 2017/745 (Class IIa), EN 455 Part 2, EN 61331-1
Quality / Environmental Standards	ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	5 pairs/box; 1 box/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## ORDERING INFORMATION

Size	6.5	7	7.5	8	8.5	9
Product Code	20873565	20873570	20873575	20873580	20873585	20873590

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 844 494 7854  
Email: insidesalesus@ansell.com  
Sandel US Tel: 866 764 3327  
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico  
Tel: +52 442 296 2050  
Email: latam.medical@ansell.com  
Ansell Brasil Ltda  
+55 (11) 3356-3100  
Email: latam.medical@ansell.com  
Ansell Colombia SAS  
Tel: +57 1 288 3247  
Email: latam.medical@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu  
Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

India  
Tel: +91 22 38124500  
Email: info.india@ansell.com  
Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com  
China  
Tel: +86(21) 3827 5005  
Email: infochina@ansell.com

**ansell.com**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2025 Ansell Limited. All Rights Reserved.